


Firmly enclosed

Reusable stone catcher

Benefit from the advantages of our stone catchers



Can be dismantled for perfect cleaning

The endoscope can be removed with the basket in situ and the handle subsequently reattached

Metal head for protection of the endoscope

Soldered for a permanent bond

Multifil




Monofil



Stainless steel

Focus on variety & efficiency.

Stone catchers are used for the endoscopic removal of gallstones and other foreign bodies. We offer our stone catchers in a range of variations.




134°C Autoclave 273°F Reusable

4-way stone catcher, **twisted**

Not suitable for lithotripsy!

Article number	Working channel	Basket diameter	Basket length in extended state	Working length	Design	PU
104.318	2.0 mm	- 10 mm	- 30 mm	200 cm	Multifil	1 unit
104.418	2.0 mm	- 15 mm	- 40 mm	200 cm		1 unit
104.518	2.0 mm	- 20 mm	- 50 mm	200 cm		1 unit
104.424	2.8 mm	- 15 mm	- 40 mm	200 cm		1 unit
104.524	2.8 mm	- 20 mm	- 50 mm	200 cm		1 unit
104.624	2.8 mm	- 25 mm	- 60 mm	200 cm		1 unit
104.724	2.8 mm	- 30 mm	- 70 mm	200 cm		1 unit
114.318	2.0 mm	- 10 mm	- 30 mm	200 cm		Monofil
114.418	2.0 mm	- 15 mm	- 40 mm	200 cm	1 unit	
114.518	2.0 mm	- 20 mm	- 50 mm	200 cm	1 unit	
114.424	2.8 mm	- 15 mm	- 40 mm	200 cm	1 unit	
114.524	2.8 mm	- 20 mm	- 50 mm	200 cm	1 unit	
114.624	2.8 mm	- 25 mm	- 60 mm	200 cm	1 unit	
114.724	2.8 mm	- 30 mm	- 70 mm	200 cm	1 unit	




134°C Autoclave 273°F

Reusable

4-way stone catcher, **roughened**

Not suitable for lithotripsy!


Article number	Working channel	Basket diameter	Basket length in extended state	Working length	Design	PU
204.318	2.0 mm	- 10 mm	- 30 mm	200 cm	Multifil	1 unit
204.418	2.0 mm	- 15 mm	- 40 mm	200 cm		1 unit
204.518	2.0 mm	- 20 mm	- 50 mm	200 cm		1 unit
204.618	2.0 mm	- 25 mm	- 60 mm	200 cm		1 unit
204.718	2.0 mm	- 30 mm	- 70 mm	200 cm		1 unit
204.424	2.8 mm	- 15 mm	- 40 mm	200 cm		1 unit
204.524	2.8 mm	- 20 mm	- 50 mm	200 cm		1 unit
204.624	2.8 mm	- 25 mm	- 60 mm	200 cm		1 unit
204.724	2.8 mm	- 30 mm	- 70 mm	200 cm		1 unit
214.318	2.0 mm	- 10 mm	- 30 mm	200 cm		Monofil
214.418	2.0 mm	- 15 mm	- 40 mm	200 cm	1 unit	
214.518	2.0 mm	- 20 mm	- 50 mm	200 cm	1 unit	
214.618	2.0 mm	- 25 mm	- 60 mm	200 cm	1 unit	
214.718	2.0 mm	- 30 mm	- 70 mm	200 cm	1 unit	
214.424	2.8 mm	- 15 mm	- 40 mm	200 cm	1 unit	
214.524	2.8 mm	- 20 mm	- 50 mm	200 cm	1 unit	
214.624	2.8 mm	- 25 mm	- 60 mm	200 cm	1 unit	
214.724	2.8 mm	- 30 mm	- 70 mm	200 cm	1 unit	



Reusable 6-way stone catcher, **twisted**

Not suitable for lithotripsy!



Article number	Working channel	Basket diameter	Basket length in extended state	Working length	Design	PU
106.424	2.8 mm	- 15 mm	- 40 mm	200 cm	Multifil	1 unit
106.524	2.8 mm	- 20 mm	- 50 mm	200 cm		1 unit
106.624	2.8 mm	- 25 mm	- 60 mm	200 cm		1 unit
106.724	2.8 mm	- 30 mm	- 70 mm	200 cm		1 unit
116.318	2.0 mm	- 10 mm	- 30 mm	200 cm	Monofil	1 unit
116.418	2.0 mm	- 15 mm	- 40 mm	200 cm		1 unit
116.518	2.0 mm	- 20 mm	- 50 mm	200 cm		1 unit
116.424	2.8 mm	- 15 mm	- 40 mm	200 cm		1 unit
116.524	2.8 mm	- 20 mm	- 50 mm	200 cm		1 unit
116.624	2.8 mm	- 25 mm	- 60 mm	200 cm		1 unit
116.724	2.8 mm	- 30 mm	- 70 mm	200 cm		1 unit



Reusable 6-way stone catcher, **roughened**

Not suitable for lithotripsy!

Article number	Working channel	Basket diameter	Basket length in extended state	Working length	Design	PU
206.424	2.8 mm	- 15 mm	- 40 mm	200 cm	Multifil	1 unit
206.524	2.8 mm	- 20 mm	- 50 mm	200 cm		1 unit
206.624	2.8 mm	- 25 mm	- 60 mm	200 cm		1 unit
206.724	2.8 mm	- 30 mm	- 70 mm	200 cm		1 unit
216.318	2.0 mm	- 10 mm	- 30 mm	200 cm		Monofil
216.418	2.0 mm	- 15 mm	- 40 mm	200 cm	1 unit	
216.518	2.0 mm	- 20 mm	- 50 mm	200 cm	1 unit	
216.424	2.8 mm	- 25 mm	- 40 mm	200 cm	1 unit	
216.524	2.8 mm	- 20 mm	- 50 mm	200 cm	1 unit	
216.624	2.8 mm	- 25 mm	- 60 mm	200 cm	1 unit	
216.724	2.8 mm	- 30 mm	- 70 mm	200 cm	1 unit	



Reusable 8-way stone catcher, **twisted**

Not suitable for lithotripsy!


Article number	Working channel	Basket diameter	Basket length in extended state	Working length	Design	PU
108.424	2,8 mm	- 15 mm	- 40 mm	200 cm	Multifil	1 unit
108.524	2,8 mm	- 20 mm	- 50 mm	200 cm		1 unit
108.624	2,8 mm	- 25 mm	- 60 mm	200 cm		1 unit
108.724	2,8 mm	- 30 mm	- 70 mm	200 cm		1 unit
118.424	2,8 mm	- 15 mm	- 40 mm	200 cm		Monofil
118.524	2,8 mm	- 20 mm	- 50 mm	200 cm	1 unit	
118.624	2,8 mm	- 25 mm	- 60 mm	200 cm	1 unit	
118.724	2,8 mm	- 30 mm	- 70 mm	200 cm	1 unit	




Reusable 8-way stone catcher, **roughened**

Not suitable for lithotripsy!

Article number	Working channel	Basket diameter	Basket length in extended state	Working length	Design	PU
208.424	2,8 mm	- 15 mm	- 40 mm	200 cm	Multifil	1 unit
208.524	2,8 mm	- 20 mm	- 50 mm	200 cm		1 unit
208.624	2,8 mm	- 25 mm	- 60 mm	200 cm		1 unit
208.724	2,8 mm	- 30 mm	- 70 mm	200 cm		1 unit
218.424	2,8 mm	- 15 mm	- 40 mm	200 cm		Monofil
218.524	2,8 mm	- 20 mm	- 50 mm	200 cm	1 unit	
218.624	2,8 mm	- 25 mm	- 60 mm	200 cm	1 unit	
218.724	2,8 mm	- 30 mm	- 70 mm	200 cm	1 unit	



More products from Maslanka



Company

Since 1974, **H. + H. Maslanka Chirurgische Instrumente GmbH** has been known for its high-quality endoscopy equipment. Thanks to our comprehensive portfolio, we offer a wide range of products and cover the entire field of endoscopy.

Our customers include medical technology companies, pharmaceutical companies, as well as hospitals and practices.

Products for gastroenterology, cardiology and urology as well as special solutions for surgery are all part of our product range.

In our product catalogue you will find a detailed overview of our products. Please do not hesitate to contact us if you have any further questions.

- Gastroenterology
- Cardiology
- Urology
- Neurology
- Customer-specific developments

Reusable Stone catcher



H. + H. Maslanka
Chirurgische Instrumente GmbH
Stockacherstraße 172
78532 Tuttlingen / Germany

Tel.: +49 7461 9607 - 0
Fax: +49 7461 - 47 32
info@maslanka.de
www.maslanka.de

07/2021 | EN

CE 0483