



Multi-use accessories
Handles & high-frequency cables




H. + H. Maslanka
Chirurgische Instrumente GmbH
Stockacherstraße 172
78532 Tuttlingen / Germany

Tel.: +49 7461 9607 - 0
Fax: +49 7461 - 47 32
info@maslanka.de
www.maslanka.de

Matching Accessories

Multi-use handles and high-frequency cables

CE 0483 09/2020 | EN



134°C Autoclave 273°F Multi-use handle

The right handle for our products.

Article number	Connection	Suitable for	PU
6202LL	Short handle with HF connector for Olympus cable	High-frequency forceps	1
6203LL	Long handle with HF connector for Olympus cable	Loops	1



Gold-plated for better contact



6403A, 4 mm
6405A, 5 mm
6400LL, 4 mm
6404LL, 5 mm



6406A, 6 mm
6402LL, 6 mm



6408A, 8 mm
6407LL, 8 mm

134°C Autoclave 273°F Multi-use high-frequency cable

Article number	Connection handle side	Device side	Working length
6400LL	for HM handle	4 mm plug, Martin, Storz, Erbe (E and T series)	300 cm
6404LL	for HM handle	5 mm plug, Erbe (ICC device and VIO devices)	
6402LL	for HM handle	6 mm plug, Olympus (120/300 device and VIO devices)	
6407LL	for HM handle	8 mm plug (Vallylab)	
6403A	Olympus	4 mm plug, Martin, Storz, Erbe (E and T series)	
6405A	Olympus	5 mm plug, Erbe (ICC device and VIO devices)	
6406A	Olympus	6 mm plug, Olympus (120/300 device and VIO devices)	
6408A	Olympus	8 mm plug, Vallylab	